

## KNEE BRACE FITTING STEP-BY-STEP

1. Open brace – wide part on top.



2. **A)** Fasten upper thigh cuff.  
**B)** Upper thigh cuff fastened.



3. **A)** Fasten lower calf cuff.  
**B)** Lower calf cuff fastened.



4. Pull and fasten upper thigh strap for comfortable compression.



## KNEE BRACE FITTING STEP-BY-STEP

5. Pull and fasten lower calf strap for comfortable compression.



6. Hinges should be anchored on both sides of the knee joint. Double check fit once more to ensure proper tension.

