



Cruiser III Wheelchair Lightweight, Dual-Axle Wheelchair

HCPCS: K0003



FIGURE A



FIGURE B



FIGURE C



FIGURE D




FIGURE E

- Built-in seat rail extensions and extendable upholstery easily adjust seat depth from 16" to 18" (Figure A - C)
- Weighs under 36 lbs (excluding front riggings)
- Carbon steel frame with silver vein finish
- Removable flip-back arms allow for easy transfer (Figure D)
- New frame style eliminates seat guides and allows for custom back inserts and accessories
- Nylon upholstery is durable, lightweight, attractive and easy-to-clean
- Composite, Mag-style wheels are lightweight and maintenance free
- 8" front casters are adjustable in three positions
- Padded armrests provide additional comfort
- Comes with swing-away footrests or elevating leg rests with tool-free adjustable length riggings (Figure E)
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Dual axle provides easy transition of seat height to hemi-level
- Comes with push-to-lock wheel locks

PRODUCTS

Description

PRODUCTS	Description
K316DDA-SF	Cruiser III 16Lt.Wt FlipBk RemArm, 1/cs
K316DDA-ELR	Cruiser III 16Lt.Wt FlipBk DeskArm, 1/cs
K316DFA-SF	Cruiser III 16Lt.Wt FlipBk FullArm, 1/cs
K316DFA-ELR	Cruiser III 16Lt.Wt FlipBk FullArm, 1/cs
K316ADDA-ELR	Cruiser III 16Lt.Wt FlipBkDskArm, 1/cs
K316ADFA-ELR	Cruiser III 16Lt.Wt FlipBkFullArm, 1/cs
K316ADFA-SF	Cruiser III 16Lt.Wt FlipBkFullArm, 1/cs
K316ADDA-SF	Cruiser III 16Lt.Wt FlipBk, 1/cs
K318DDA-SF	Cruiser III 18Lt.Wt FlipBk SwgFt, 1/cs
K318DDA-ELR	Cruiser III 18Lt.Wt FlipBk ElvLgRst, 1/cs
K318DFA-ELR	Cruiser III 18Lt.Wt FlipBk FullArm, 1/cs
K318DFA-SF	Cruiser III 18Lt.Wt FlipBk, 1/cs
K318ADDA-SF	Cruiser III 18 Wt Flip Back, 1/cs
K318ADFA-SF	Cruiser III 18Lt.Wt FlipBkFullArm, 1/cs
K318ADFA-ELR	Cruiser III 18Lt.Wt FlipBkFullArm, 1/cs
K318ADDA-ELR	Cruiser III 18Lt.Wt FlipBk DeskArm, 1/cs
K320DDA-SF	CruiserIII 20Lt.Wt FlipBk RemDskArm, 1/cs
K320DDA-ELR	Cruiser III 20Lt.Wt FlipBk DskArm, 1/cs
K320DFA-SF	Cruiser III 20Lt.Wt FlipBk FullArm, 1/cs
K320DFA-ELR	Cruiser III 20Lt.Wt FlipBk FullArm, 1/cs
K320ADDA-SF	Cruiser III 20Lt.Wt FlipBk RemArm, 1/cs
K320ADDA-ELR	Cruiser III 20Lt.Wt FlipBk, 1/cs
K320ADFA-ELR	Cruiser III 20Lt.Wt FlipBkFullArm, 1/cs

PRODUCTS	Description
K320ADFA-SF	Cruiser III 20Lt.Wt FlipBkFullArm, 1/cs
SPECIFICATIONS	Description
Width Closed	12"
Height	35"
Seat--To-Floor Height	17.5" - 19.5"
Back Height	16"
Seat Depth	16" - 18"
Length with Riggings	42"
Product Weight (with 16" Seat)	35 lbs.
Product Weight (with 18" Seat)	36 lbs.
Product Weight (with 20" Seat)	38 lbs.
Weight Limit (with 20" seat)	350 lbs.
Weight Limit (with 16" and 18" Seat)	300 lbs.
Carton Shipping Weight	43 lbs. or less
Warranty	Limited Lifetime
ACCESSORIES	Description
STDS801GT	Wheel Lock Extension, 1/ea
STDS801	Wheel Lock Extension, 1/ea
STDS1002	Standard 8" Caster, 1/ea
STDS806	Anti-Fold Bar, 1/ea
STDS819	Anti-Tippers With Wheels, 1/pr
STDS831	Heel Loops, 1/bx
WASR	Limb Support, Right, 1/ea
STDS834	Overhead Anti-Theft Device Single Pole, 1/ea
STDS3J24SF	Plastic Foot plate, 1/pr
LK3JELR	ELR Set, 1/pr
STDS850	Seat Belts, Auto-Clasp Type Closure, 1/ea
STDS851	Seat Belts, Hook and Loop Fastener Type Closure, 1/ea
STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric, 1/ea
STDS856	Seat Belts, Hook and Loop Fastener Type Closure, Bariatric, 1/ea
STDS820	Telescoping IV Pole Attachment, 1/ea
STDS860	Universal Forward Stabilizer, 1/pr
STDS1034	Wheelchair Cane/Crutch Holder, 1/ea
STDSSHK-K3	Super Hemi-Kit, 1/ea
STDSK234-PTL	Pull to Lock Brakes(left), 1/ea
STDSK2324-PTR	Pull to Lock Brakes(right), 1/ea
PARTS	Description
STDS801	Wheel Lock Extension, 1/ea
STDS1002	Standard 8" Caster, 1/ea
STDS806	Anti-Fold Bar, 1/ea
STDS819	Anti-Tippers With Wheels, 1/pr
STDS835	1/ea
STDS831	Heel Loops, 1/bx
WASR	Limb Support, Right, 1/ea
STDS823	Overhead Anti-Theft Device, 1/ea
STDS834	Overhead Anti-Theft Device Single Pole, 1/ea
STDS3J24SF	Plastic Foot plate, 1/pr
LK3JELR	ELR Set, 1/pr
STDS850	Seat Belts, Auto-Clasp Type Closure, 1/ea
STDS851	Seat Belts, Hook and Loop Fastener Type Closure, 1/ea
STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric, 1/ea
STDS856	Seat Belts, Hook and Loop Fastener Type Closure, Bariatric, 1/ea